# BIOFLOAT™

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# BIOFLOAT™

## Cell culture surface for spheroid culture

In many areas of biomedical research, *in vitro* models are indispensable. The most conventional form is the two-dimensional cell culture. However, discrepancies often occur in the transfer of the results to an entire organism. The aim of three-dimensional cell culture is therefore to close this gap between the *in vitro* and *in vivo* situation. Moreover, the use of three-dimensional cell models can help to overcome hurdles in the transferability of experimental results to different species and thereby reduce animal experiments.

Spheroid culture, as a simple variant of 3D cell culture, offers a reliable *in vitro* model for biomedical research. Given the three-dimensional structure, the establishment of pronounced cell-cell and cell-matrix interactions is facilitated – realistically imitating the *in vivo* microenvironment.

Spheroid cultures are a widely used culture system in drug development, toxicological studies as well as cancer and stem cell research, thereby improving efficiency and reliability in these phases of preclinical research.

Our SARSTEDT BIOFLOAT™ cell culture surface offers fast and highly reproducible generation of perfectly round-shaped spheroids. Formation of a single, symmetrical spheroid per well of our SARSTEDT BIOFLOAT™ 96-well plate qualifies BIOFLOAT™ as an excellent choice for high throughput analyses. Furthermore, the SARSTEDT BIOFLOAT™ 96-well plate satisfies the standardized dimensions (ANSI/SLAS standards) and is therefore suitable for automated processes.

The following protocol describes exemplary working steps to support the establishment of spheroid culture. Exact growth medium volume and cell seeding density depend on cell line/type and the type of downstream analyses. Therefore, careful optimization of growth and analysis conditions is recommended.

### Protocol for spheroid generation and culture


Cell seeding density depends on cell line/type, duration of spheroid cultivation and desired spheroid dimensions at the time of analysis. Taking these parameters into account, the recommended working volume per well of a 96-well plate ranges from 100 to 200  $\mu$ l.

1. Detach your cells from the culturing surface using your regular protocol (trypsinization, EDTA etc.). To seed cells for spheroid culture, a single cell suspension is recommended. If necessary, cell aggregates may be removed by straining the cell suspension through a cell strainer (83.3945. XXX).
2. Adjust the cell concentration according to the desired seeding density. The recommended density ranges from 1,000 to 10,000 cells per well (see Fig. 1).
3. Distribute the cell suspension into the wells of the BIOFLOAT™ 96-well plate (recommended working volume: 100-200  $\mu$ l). A centrifugation step is not recommended. Note that centrifugation may even negatively affect spheroid formation using BIOFLOAT™.
4. Incubate the plate in a humidified incubator (37°C and 5% CO<sub>2</sub> for most mammalian cells). For most cell types, spheroid formation in the SARSTEDT BIOFLOAT™ plate is accomplished within 2 to 24 hours. Carefully monitor spheroid formation and growth under the microscope on a daily basis.
5. Depending on duration of spheroid cultivation and cell line/type, continuous nutrient supply needs to be ensured by (partial) growth medium exchange. Removal or addition of growth medium should be performed at the side of the well to prevent spheroid damage. Furthermore, it is recommended not to remove the complete volume of growth medium, but leave a residual 10 to 20  $\mu$ l in the well. The anti-adhesive properties of the BIOFLOAT™ surface are not damaged by medium exchanges or accidental contact of the surface with the pipet tip.


## Examples of spheroid cultures in SARSTEDT BIOFLOAT™ cell culture plates


**Fig. 1:** 3T3 fibroblasts were seeded onto BIOFLOAT™ cell culture plates at different densities. As a control, a non-coated plate was used. After three days, successful spheroid formation in the BIOFLOAT™-coated plate could be confirmed by microscopic imaging, whereas the fibroblasts adhered to the non-coated control plate and failed to form spheroids. Furthermore, it could be demonstrated that keeping culturing time constant, the spheroid size strongly depends on the initial seeding density. Scale bar = 400 µm.

**Fig. 2:** Primary hepatocytes (*C. lupus familiaris*) were seeded onto BIOFLOAT™ 96-well plates (2,500 cells per well) and microscopically monitored for spheroid formation over the course of eight days. After spheroid formation, 50 µl of growth medium was exchanged every 48 - 72 h. Scale bar = 200 µm.

**Fig. 3:** Immortalized human astrocytes were seeded onto BIOFLOAT™ 96-well plates at a density of 10,000 cells/well. Subsequently, cell viability was determined over several days using fluorescence staining. Viable cells were labelled with calcein (green) and dead cells with propidium iodide (red). Scale bar = 200 µm.

**Fig. 4:** BIOFLOAT™-induced co-culture spheroids containing osteoblasts (red) and HUVEC cells (green). Scale bar = 250 µm.